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腸胃道內視鏡處理的高層次消毒的標準是由 SGNA (Society of Gastroenterology Nurses and Associates)、ASGE (American Society for Gastrointestinal Endoscopy)、ACG (The American College of Gastroenterology)、AGA (The American Gastroenterological Association)、APIC (The Association for Professionals in Infection Control and Epidemiology) 及 ASTM (American Society for Testing and Materials) 所認定。美國 CDC (The Centers for Disease Control and Prevention,) 及 JCAHO (The Joint Commission on Accreditation of Health Care Organizations) 認定高層次消毒適用於腸胃道內視鏡。只有在手術室使用內視鏡時需在無菌環境。高層次消毒可破壞所有具活動力微生物，但不包括細菌孢子。使用高層次消毒劑須依據操作指南，高層次消毒劑可重複使用，必須測試其仍保有最低有效濃度 (Minimum

effective concentration, MEC)。

(七) 高層次消毒劑 (**High-level disinfection agents**)

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- disinfection techniques eliminate the risk of endoscopic transmission of *Helicobacter pylori*. *Am J Gastroenterol* 1995; 90:227-32.
5. Spach DH, Silverstein FE, Stamm WE. Transmission of infection by gastrointestinal endoscopy and bronchoscopy. *Ann Int Med* 1993;118:117-28.
  6. Tandon RK. Endoscopic disinfection: practices and recommendations. *J Gastroenterol Hepatol* 1991;6:37-9.
  7. Rutala WA. APIC guideline for selection and use of disinfectants. *Am J Infect Control* 1990;18:99-117.
  8. Axon ATR, Cotton PB. Endoscopy and infection. *Gut* 1983;24:1064-6.
  9. Kimmey MB, Burnett DA, Carr-Locke DL, et al. Transmission of infection by gastrointestinal endoscopy. *Gastrointest Endosc* 1993; 39:885.
  10. Bond WW. Overview of infection control problems: principles in gastrointestinal endoscopy. *Gastrointest Endosc Clin N Am* 2000; 10:199.
  11. Nelson DB. Infectious disease complications of GI

- endoscopy: part II, exogenous infections. *Gastrointest Endosc* 2003; 57:695.
12. Wallace CG, Agee PM, Demicco DD. Liquid chemical sterilization using peracetic acid. An alternative approach to endoscope processing. *ASAIIO J* 1995; 41:151.
  13. Kovacs BJ, Chen YK, Kettering JD, et al. High-level disinfection of gastrointestinal endoscopes: are current guidelines adequate? *Am J Gastroenterol* 1999; 94:1546.
  14. Foliente RL, Kovacs BJ, Aprecio RM, et al. Efficacy of high-level disinfectants for reprocessing GI endoscopes in simulated-use testing. *Gastrointest Endosc* 2001; 53:456.
  15. Epstein L, Hunter JC, Arwady MA, et al. New Delhi metallo- $\beta$ -lactamase-producing carbapenem-resistant *Escherichia coli* associated with exposure to duodenoscopes. *JAMA* 2014; 312:1447.
  16. Rutala WA, Weber DJ. Gastrointestinal endoscopes: a need to shift from disinfection to sterilization? *JAMA* 2014; 312:1405.

17. ASGE Quality Assurance In Endoscopy Committee, Petersen BT, Chennat J, et al. Multisociety guideline on reprocessing flexible gastrointestinal endoscopes: 2011. *Gastrointest Endosc* 2011; 73:1075.
18. Nelson DB, Barkun AN, Block KP, et al. Technology status evaluation report. Transmission of infection by gastrointestinal endoscopy. May 2001. *Gastrointest Endosc* 2001; 54:824.
19. American Society for Gastrointestinal Endoscopy, Petersen BT, Adler DG, et al. Automated endoscope reprocessors. *Gastrointest Endosc* 2009; 69:771.
20. Jung M, Beilenhoff U, Pietsch M, et al. Standardized reprocessing of reusable colonoscopy biopsy forceps is effective: results of a German multicenter study. *Endoscopy* 2003; 35:197.
21. Furman PJ. Third-party reprocessing of endoscopic accessories. *Gastrointest Endosc Clin N Am* 2000; 10:385.
22. Systchenko R, Marchetti B, Canard JN, et al. Guidelines of the French Society of Digestive Endoscopy:

- recommendations for setting up cleaning and disinfection procedures in gastrointestinal endoscopy. *Endoscopy* 2000; 32:807.
23. Society of Gastroenterology Nurses and Associates, Inc. Standards of infection control in reprocessing of flexible gastrointestinal endoscopes. *Gastroenterol Nurs* 2000; 23:172.
24. American Society for Gastrointestinal Endoscopy, Petersen BT, Adler DG, et al. Automated endoscope reprocessors. *Gastrointest Endosc* 2009; 69:771.
25. Humphreys H, McGrath H, McCormick PA, Walsh C. Quality of final rinse water used in washer-disinfectors for endoscopes. *J Hosp Infect* 2002; 51:151.
26. Ahuja V, Tandon RK: Survey of gastrointestinal endoscope disinfection and accessory reprocessing practices in the Asia-Pacific region. *J Gastroenterol Hepatol* 2000;15:78-81.
27. Alfa MJ, Olson N, Degagné P, Simner PJ. Development and validation of rapid use scope test strips to determine

the efficacy of manual cleaning for flexible endoscope channels. *Am J Infect Control* 2012; 40:860.

28. Bommarito M, Thornhill GA, Morse DJ. A Multi-site Field Study Evaluating the Effectiveness of Manual Cleaning of Flexible Endoscopes with an ATP Detection System. *Am J Infect Control* 2013; 41:S24.